



Alma Medical

www.alma-medical.eu

MEDICAL PRODUCTS CATALOG



DRESSING MATERIALS



ABSORBENT GAUZE

5 cm x 5 cm	7.5 cm x 15 cm
7.5 cm x 7.5 cm	10 cm x 10 cm
8 cm x 18 cm	10 cm x 20 cm

Used as a wound covering or as a mechanical barrier for absorption or suppression at surgical sites and wounds.



HYDROPHILIC GAUZE

Width: 90 cm
Length: 100 m

Hydrophilic gauze is used as a mechanical barrier, as a wound covering, and for absorption/suppression during surgical procedures.

Composition: 100% cotton yarn



GAUZE SWABS

Used as a wound covering or as a mechanical barrier for absorption or suppression at surgical sites and wounds.



STERILE COMPRESSES

30 cm x 30 cm	40 cm x 40 cm
45 cm x 45 cm	

Used in abdominal surgical procedures.



DRESSING MATERIALS



PLASTER OF PARIS BANDAGE

Width: 5 cm	Width: 7 cm
Width: 10 cm	Width: 15 cm
Width: 20 cm	

Made from 100% cotton and orthopedic plaster.
Fast setting: 2-8 minutes.
Easy to shape and adjust.



HYDROPHILIC BANDAGE

Width: 5 cm	Width: 7 cm
Width: 10 cm	Width: 15 cm
Width: 20 cm	

Made with firm knitted edges that do not loosen, reducing the risk of infection. Provides secure, hygienic and comfortable application during every medical procedure.



ELASTIC BANDAGE

6 cm	8 cm	10 cm
12 cm	15 cm	20 cm
2 m	2.5 m	3 m

Used for fixation of minor injuries requiring compression, to stabilize muscles during recovery from muscle injuries, and to reduce strain during movement.



CAST PADDING

6 cm	8 cm	10 cm
12 cm	15 cm	20 cm
2 m	2.5 m	3 m

Designed to prevent irritation where the plaster bandage is wrapped. Made from medical-grade synthetic fibers. Helps prevent and reduce wound friction and irritation.



DRESSING MATERIALS



HYDROPHILIC COTTON PAD

Medical material made from gauze produced with hydrophilic, highly absorbent cotton yarn.

Used in surgical procedures, for open wounds and abdominal operations, and as a post-partum dressing material.



HYDROPHILIC COTTON WOOL

Made from 100% cotton

Manufactured in accordance with TS 4786

Non-sterile

Pack of 1,000 g



TRANSPARENT TAPE

1.25 cm x 5 m, boxed
2.5 cm x 5 m, boxed

Recommended for dialysis patients and whenever monitoring the skin is important, while providing protection and comfort.

TRANSPARENT STRUCTURE - Protects the wound from water and other contamination.

HYPOALLERGENIC - Does not cause skin irritation.

BREATHABLE - Allows the skin to breathe.

TRANSPARENT - Keeps the wound and skin visible.

PRACTICAL - Easy to tear lengthwise and crosswise.

RADIOTRANSSPARENT - X-rays can be taken without removing the tape.

VARIOUS SIZES - Available in different sizes and packs.



DRESSING MATERIALS



SURGICAL PAPER TAPE

2.5 cm x 5 m, boxed
5 cm x 5 m, boxed

The microporous structure allows free airflow and natural moisture transfer, keeping the skin healthy. A dispenser is also available for faster and easier application.

The tape is radiotransparent, so it does not need to be removed for X-rays. It is available in several sizes and pack formats for different needs.



SURGICAL FABRIC TAPE

1.5 cm x 5 m, boxed
2.5 cm x 5 m, boxed
5 cm x 5 m, boxed

Made from 100% woven cotton fabric in a neutral color, with a shelf life of 5 years.

The adhesive strength is at least 200 gf/cm. Manufactured in accordance with TS 3957.



DRESSING MATERIALS



NONWOVEN WOUND DRESSINGS

9 cm x 10 cm
25 x 1

A special absorbent pad in the center absorbs blood and other wound fluids and protects the wound from external factors.

Radiolucent, allowing X-rays without removal.

Sterilized by gamma irradiation.

Remains sterile until the package is opened or damaged.



EYE PATCHES

Adults: 5.8 x 8.3 cm (5 x 1)
Children: 5.7 x 7.3 cm (5 x 1)

Intended for children undergoing treatment for amblyopia (lazy eye) or strabismus (crossed eyes).

The black backing covers the healthy eye to strengthen the treated eye.



SPECIMEN COLLECTION MATERIALS



STERILE URINE CONTAINER

Special design with a hermetic seal.

60 ml 100 ml



STOOL CONTAINER

Polypropylene tube with a volume of 15 ml.



URINE CONTAINER (3,000 ml)

Stores urine for up to 24 hours. Special cap for transferring the specimen.



ESMARCH TOURNIQUET

Easy to apply and remove, with a tightly woven flexible strap.



3-PLY MASK WITH EAR LOOPS

Pack of 50 pcs.



WHITE KN95 MASK

Made from high-quality materials. 95% protection against toxic dust. Adapts to every face shape.



PROTECTIVE EQUIPMENT



DISPOSABLE COVERALL

Protective coverall for single use.



DISPOSABLE GOWN

Made from nonwoven material; comfortable and practical. Provides hygienic protection during medical and sanitary procedures.



DISPOSABLE APRON

Hypoallergenic, hygienic and fluid-resistant aprons that tie easily around the waist.

Available in various colors.



PVC SHOE COVERS (1,000 PCS.)

Non-sterile. Made from polypropylene.



NONWOVEN SHOE COVERS (100 PCS.)

Non-sterile. Made from nonwoven fabric.



WHITE/BLUE CAPS (100 PCS.)

Made from nonwoven fabric.



PVC SLEEVE COVERS

Waterproof, with elastic cuffs.



PLASTERS



ADHESIVE PLASTER 19 x 72 mm (10/1)

Adhesive plaster for minor cuts, suitable for any part of the body. The zinc oxide adhesive retains its adhesion for a long time without irritation.



ADHESIVE PLASTER PLUS 19 x 72 mm (10/1)

With a special pad, practical carrier and unique design, this product provides lasting adhesion without skin irritation.



TRADITIONAL PLASTERS (20/1)

For those who prefer classic first-aid plasters for minor wounds; a classic plaster with a special pad.



WATERPROOF PLASTERS 19 x 72 mm (10/1)

Ideal for daily use; protects wounds during contact with water.



UNIVERSAL MIXED PLASTERS (20/1)

19 x 72 mm	16 x 55 mm
25 x 72 mm	Ø25.5 mm



PLASTERS FOR SENSITIVE SKIN (20/1)

Specially designed for sensitive skin, with nonwoven material and a safe adhesive; ideal for older adults and children.



PLASTERS



PLASTERS FOR LARGER CUTS 26 x 72 mm (10/1)

The special pad helps stop bleeding from larger cuts that do not require surgical treatment.



ALOE VERA PLASTERS (20/1)

Protective plasters for minor wounds and scratches. Suitable for general use on small cuts; helps reduce skin irritation and provides a cooling sensation.



CAPSICUM PLASTER 10 x 18 cm

Perforated Cansinplast heat plaster for rheumatism, sciatica, lumbago, muscle pain, and back or neck spasms. Relaxes the application area by providing warmth.



BACK PAIN PLASTER 10 x 13 cm

Warms the application area for 10-12 hours at approximately 45-50°C. Apply only to clothing, not directly to the body.



HOT & COLD GEL - 13 x 26 cm / 13 x 13 cm

HOT COMPRESS: To increase blood circulation; for rheumatism, muscle pain, muscle stiffness, back pain and neck muscle spasms.

COLD COMPRESS: For sports injuries or bruising; fractures, swelling and muscle tears; headache and migraine; nosebleeds; toothache; and first aid for burns and insect bites.



PLASTERS



ANTI-SNORING NASAL STRIPS 67 x 20 mm (10/1)

For allergic congestion, swelling of the nasal septum, a blocked nose caused by a cold, or difficult breathing.



MOSQUITO-REPELLENT PATCHES 15 x 24 mm (18/1)

The fragrance of eucalyptus, citronella and geranium essential oils keeps mosquitoes away. Suitable for children over 36 months. Effective for 8 hours; apply only to clothing.



ACNE PATCHES Ø15 mm (20/1)

Helps prevent infection of non-inflammatory acne and absorbs secretions from inflamed spots. Waterproof, breathable and hypoallergenic. Apply to clean, dry skin before makeup.



BLACKHEAD REMOVAL STRIPS (4/1)

Apply to areas with blackheads and leave for 15 minutes. Remove the strip; blackheads and impurities remain on it, cleansing the application area.



BREATHING RELIEF PATCHES (3/1)

Made with natural eucalyptus, mint and pine extracts. Apply to clothing; the pleasant aroma soothes, supports easier breathing and helps clear the nose.



TRANSPARENT ADHESIVE PLASTER

Ultra-thin and discreet on the skin. Ideal for minor injuries.



NEEDLES & CANNULAS



INJECTION NEEDLES

■	18G 38 mm - 18G x 1½" (1.20 x 38 mm)
■	20G 38 mm - 20G x 1½" (0.90 x 38 mm)
■	21G 38 mm - 21G x 1½" (0.80 x 38 mm)
■	22G 32 mm - 22G x 1¼" (0.70 x 32 mm)
■	23G 32 mm - 23G x 1¼" (0.60 x 32 mm)
■	26G 13 mm - 26G x ½" (0.45 x 13 mm)



TWO-PIECE SYRINGE

5 ml x 250	2 ml x 350
10 ml x 150	20 ml x 100

THREE-PIECE SYRINGE

1 ml x 100	2 ml x 350	5 ml x 250
10 ml x 200		



INFUSION SET

Sterile, single-use and non-pyrogenic;
Flow rate: 20 drops / 1 ml;
Roller clamp for adjusting the drip rate;
Closed air channel;
For serum in glass bottles and plastic bags;
Filter chamber;
Length: 150 cm
Needle: 21G




INTRAVENOUS CANNULA

16G (gray)	18G (green)	20G (pink)
22G (blue)	24G (yellow)	26G (purple)



MINIVEN (BABY SYSTEM)

21 (green)	22 (black)	23 (blue)
24 (purple)		25 (orange)



SUTURES



SUTURES FOR WOUND CLOSURE

Used during surgery or for injuries to provide strong, secure wound closure, prevent infection and reduce the risk of complications.

ABSORBABLE SUTURE

This product is widely used in neurosurgery and cardiovascular surgery, where high precision and reliability are required during treatment and procedures.

It also helps reduce the risk of infection and prevent inflammatory reactions, supporting safer and more effective patient recovery.



NON-ABSORBABLE SUTURE

Specially designed for microscopic and delicate surgical procedures requiring maximum precision and stability.

Its properties provide high resistance to infections that may occur in or around the surgical area, reducing complications and enabling faster recovery.


High strength prevents elongation during use.



CATHETERS



FOLEY CATHETER LATEX (ALL SIZES)



Made from natural latex with a silicone coating.
Rubber valve for different requirements.

FOLEY CATHETER SILICONE (ALL SIZES)



Made from 100% medical-grade silicone.

GLOVES

LATEX GLOVES

Powdered S, M, L x 100

Powder-free S, M, L x 100



NITRILE GLOVES

Powder-free S, M, L x 100

Powder-free and latex-free,
non-allergenic, and tested for use with
chemotherapy drugs.

Textured fingertips for improved
precision; single-use.
Pack of 100, AQL 1.5.



PERSONAL HYGIENE PRODUCTS



COTTON PADS

Available in packs of 50, 70 and 100. 100% cotton.



COTTON SWABS



ZIG-ZAG COTTON WOOL

Available in packs of 50 g, 100 g and 200 g.



WET WIPES

Available in packs of 40, 70 and 100.



MEASURING INSTRUMENTS



INHALER



For the treatment of asthma, allergies and other respiratory conditions.

PULSE OXIMETER



Dual OLED digital display
SpO2 and heart-rate monitoring
4-direction, 6-mode display

DIGITAL WRIST BLOOD PRESSURE MONITOR



Precise digital measurement of blood pressure and pulse.

DIGITAL BLOOD PRESSURE MONITOR



Easy-to-use, precise device for accurate blood pressure and pulse measurement.

DIGITAL THERMOMETER



Precise measurement; automatic shut-off.
LCD display
Result in 60 seconds.

NON-CONTACT THERMOMETER



Fast, reliable body and surface temperature measurement. Backlit LCD display.

STETHOSCOPE



Precise acoustics for clear heart and lung sounds.

OTHER PRODUCTS



MALE URINAL



An aid that allows male patients to urinate from bed. Simple and easy to use.

FEMALE PVC URINAL 1,000 ml



Practical, hygienic urine collector for women. Made from durable PVC.

PLASTIC FEMALE BEDPAN



Used as an aid for urination in bed. Suitable for bedridden patients who cannot get up.

HOT-WATER BOTTLES / WITH COVER



Comfortable hot-water bottle with cover for heat or cold therapy.

OTHER PRODUCTS




WOODEN TONGUE DEPRESSORS x 100



150 mm x 18 mm

PERFORATED BED SHEET COVER



Laminated two-layer cover

45 x 50
50 x 50
60 x 50
70 x 50

ABSORBENT UNDERPAD



60 x 90 cm - 30 g

60 x 90 cm - 34 g

60 x 90 cm - 44 g

FIRST AID CABINET

CONTENTS:

Surgical fabric tape 1.25 cm x 5 m - 2 pcs.
Sterile gauze 7.5 x 7.5 cm - 2 packs (10 pcs.)
Hydrophilic bandage 5 cm x 4 m - 2 pcs.
Hydrophilic bandage 7 cm x 4 m - 2 pcs.
Hydrophilic bandage 15 cm x 4 m - 2 pcs.
Conforming elastic bandage 6 cm x 4 m - 2 pcs.
Adhesive plasters 19 x 72 mm - 2 packs (20 pcs.)
Mixed plasters - 1 pack (20 pcs.)
Zig-zag cotton wool - 50 g
Sterile postoperative dressing 9 x 10 cm - 2 pcs.



MEDICAL WASTE CONTAINERS



AND OTHER PLASTIC CONTAINERS

MODULAR BIN FOR GLOVES AND MASKS



(With wheels / without wheels) Bin and lid.

MEDICAL WASTE CONTAINER

30 L - 390 x 390 x 830 mm
40 L - 390 x 390 x 830 mm
60 L - 390 x 390 x 830 mm



Made from virgin material with cut resistance.

TROLLEY FOR TRANSPORTING HAZARDOUS WASTE



TATA-series trolley designed for the safe transport of hazardous-waste containers.

MOBILE MEDICAL WASTE SET



Container size: 390 x 390 x 1,030 mm
Bin size: 30 L, 390 x 390 x 830 mm

BLOOD SAMPLE SET



Contains 68 or 180/200 specimen slots. Designed for built-in installation.

MEDICAL WASTE CONTAINER



0.7 L container, 570 x 450 x 490 mm
1.2 L container, 570 x 450 x 490 mm
1.3 L container, 570 x 450 x 540 mm
1.5 L container, 570 x 450 x 490 mm



Alma Medical

St. Todosija Paunov no.8, 2300

Kochani North Macedonia

+389 33 271 666

info@alma-medical.eu

www.alma-medical.eu